

# Hydrogels

[DOWNLOAD HERE](#)

1;Preface;5 2;Table of Contents ;6 3;List of Contributors;8 4;Hydrogels and Tissue Engineering;11 4.1;1 Introduction;11 4.2;2 Materials;12 4.2.1;2.1 Synthetic materials;13 4.2.2;2.2 Naturally derived materials;13 4.3;3 Scaffold design;14 4.4;4 Degradation;15 4.5;5 Diffusion;16 4.6;6 Biological properties;16 4.7;References;18 5;Structure-Property Relationships in Hydrogels;19 5.1;1 Hydrogel classification and basic structure;19 5.1.1;1.1 Structural parameters;21 5.2;2 Hydrogel mechanical properties;21 5.2.1;2.1 Dynamical mechanical analysis;21 5.2.2;2.2 Hydrogel viscoelastic properties;23 5.2.3;2.3 Stress-strain behavior;24 5.3;3 Hydrogel swelling properties;25 5.4;References;28 6;Water and Surfaces: a Linkage Unexpectedly Profound;31 6.1;References;33 7;Polysaccharide Based Hydrogels for Biomedical Applications ;35 7.1;1 Hydrogels for soft tissue regeneration;36 7.1.1;1.1 Hydrogel for cartilage regeneration;36 7.1.2;1.2 Amidic Carboxymethylcellulose hydrogel;36 7.1.2.1;1.2.1 Hydrogels for nucleus pulposus replacement;43 7.1.3;1.3 Amidic alginate hydrogel;43 7.1.4;1.4 Hydrogels for delivery of biologically active substances;48 7.1.4.1;1.4.1 Hydrogels as vehicles for bioactive substances for regeneration of bone tissue ;48 7.1.4.2;1.4.2 Guar gum based Hydrogel enriched with HA;48 7.2;2 Conclusions;50 7.3;References;50 8;Hydrogels for Healing;52 8.1;1 Introduction;52 8.2;2 Sphere templated polymers;54 8.3;3 A biodegradable pHEMA hydrogel;55 8.4;4 Degradation studies;56 8.5;5 Conclusions;59 8.6;References;59 9;Stereocomplexed PEG-PLA Hydrogels;61 9.1;1 Introduction;61 9.2;2 Hydrogels;62 9.2.1;2.1 General requirements;62 9.2.2;2.2 Synthetic polymers used for hydrogel preparation;62 9.2.3;2.3 In situ forming hydrogels;63 9.2.4;2.4 Physical crosslinking;63 9.2.5;2.5 Thermosensitive hydrogels;64 9.2.6;2.6 Stereocomplexation;64 9.2.7;2.7 Chemical crosslinking by photopolymerization;65 9.2.8;2.8 Synthesis;65 9.2.9;2.9 Solubility;67 9.2.10;2.10 Gelation behavior;67 9.2.11;2.11 Rheology;69 9.2.12;2.12 In situ monitoring of mechanical properties during photopolymerization;70 9.3;3 Conclusions;70 9.4;References;71 10;Hybrid Hydrogels Based on Poly (vinylalcohol) Chitosan Blends and Relevant CNT Composites ;74 10.1;1 Introduction;74 10.2;2 Experimental procedures;76 10.2.1;2.1 Materials;76 10.2.2;2.2 Processing;76 10.2.3;2.3 Characterization;77 10.3;3 Results and discussion;78 10.4;4 Conclusions;83 10.5;References;85

11; Poloxamine Hydrogels: from low Cell Adhesion Substrates to Matrices with Improved Cytocompatibility for Tissue Engineering Applications ;86 11.1;1 Introduction;86 11.2;2 Poloxamers (Pluronicf) ;87 11.3;3 Poloxamines (Tetronic) ;88 11.3.1;3.1 Poloxamine/collagen hybrids;88 11.3.1.1;3.1.1 Design of semi-interpenetrating networks and their characterization;88 11.3.1.2;3.1.2 Cell embedding;90 11.3.1.3;3.1.3 Cell attachment on poloxamine/collagen hydrogels;90 11.3.2;3.2 Positively-charged matrices;93 11.3.2.1;3.2.1 Copolymerization of reactive poloxamine with quaternary ammonium methacrylates ;93 11.3.2.2;3.2.2 Methylation of poloxamine central block (quaternized poloxamine);94 11.3.3;3.3 The combined approach;95 11.3.4;3.4 The degradable strategy;96 11.4;4 Summary and perspectives;99 11.5;References;100 12; Biohydrogels for the In Vitro Re-construction and In Situ Regeneration of Human Skin ;104 12.1;1 Biohydrogels as scaffolds for in vitro produced skin equivalents;104 12.2;2 Biohydrogels for in situ human skin regeneration;107 12.3;3 Biohydrogels as carriers of topical drugs facilitating wound healing and attenuating cutaneous inflammation ;109 12.4;4 Anti-oxidant properties of biohydrogels loaded with plant polyphenols ;110 12.5;5 Biohydrogels for drug delivery through the cutaneous barrier;112 12.6;6 Perspectives of biohydrogels for dermatological and cosmetological applications ;113 12.7;References;113 13; Chitosan-Based Beas EAN/ISBN : 9788847011045 Publisher(s): Springer, Berlin, Springer, Milan Format: ePub/PDF Author(s): Barbucci, Rolando

[DOWNLOAD HERE](#)

Similar manuals:

[Hydrogels](#)